


CLIP ANEURISM YASARGIL STANDARD




*pentru inchidere permanenta



Cod	244-400-07	244-400-09	244-400-11	244-400-15	244-400-17	244-400-20
	244-400-07T	244-400-09T	244-400-11T	244-400-15T	244-400-17T	244-400-20T
Deschidere maxima	6,2 mm	7,0 mm	7,8 mm	9,2 mm	10,6 mm	11,4 mm
Forta inchidere	150 gf	180 gf	180 gf	200 gf	200 gf	200 gf

Titan


Otel pentru implant



Cod	244-400-25	244-401-07	244-401-09	244-401-11	244-401-15	244-401-20
	244-400-25T	244-401-07T	244-401-09T	244-401-11T	244-401-15T	244-401-20T
Deschidere maxima	13,3 mm	6,0 mm	6,8 mm	7,5 mm	8,7 mm	11,0 mm
Forta inchidere	200 gf	150 gf	180 gf	180 gf	200 gf	200 gf

Titan

Otel pt. implant




Cod	244-403-07	244-403-09	244-403-11	244-403-15	244-403-17
	244-403-07T	244-403-09T	244-403-11T	244-403-15T	244-403-17T
Deschidere maxima	5,5 mm	6,1 mm	6,5 mm	7,4 mm	9,9 mm
Forta inchidere	150 gf	180 gf	180 gf	200 gf	200 gf

Titan


Otel pt. implant

*pentru inchidere permanenta




Cod	244-405-07	244-405-09	244-405-11	244-410-07	244-410-09	244-410-12
	244-405-07T	244-405-09T	244-405-11T	244-410-07T	244-410-09T	244-410-12T
Deschidere maxima	5,4 mm	5,8 mm	6,2 mm	7,9 mm	8,7 mm	10,2 mm
Forta inchidere	180 gf	180 gf	180 gf	200 gf	200 gf	200 gf

Titan Otel pt. implant



Cod	244-410-17	244-410-20	244-411-05	244-411-07	244-411-09	244-411-11
	244-410-17T	244-410-20T	244-411-05T	244-411-07T	244-411-09T	244-411-11T
Deschidere maxima	11,9 mm	13,0 mm	6,2 mm	6,5 mm	7,8 mm	8,4 mm
Forta inchidere	180 gf	180 gf	200 gf	200 gf	200 gf	200 gf

Titan Otel pt. implant




Cod	244-413-05	244-413-07	244-413-10	244-415-08	244-415-10	244-417-05
	244-413-05T	244-413-07T	244-413-10T	244-415-08T	244-415-10T	244-417-05T
Deschidere maxima	5,6 mm	5,6 mm	5,6 mm	5,4 mm	5,6 mm	7,2 mm
Forta inchidere	200 gf	200 gf	200 gf	200 gf	200 gf	200 gf

Titan Otel pt. implant



CLIP ANEURISM YASARGIL STANDARD



*pentru inchidere permanenta



Cod	244-419-09	244-419-11	244-421-07	244-421-11	244-423-09	244-423-11
	244-419-09T	244-419-11T	244-421-07T	244-421-11T	244-423-09T	244-423-11T
Deschidere maxima	7,2 mm	8,2 mm	5,7 mm	8,0 mm	6,2 mm	6,8 mm
Fora inchidere	200 gf	200 gf	200 gf	200 gf	200 gf	200 gf


 Titan  Otel pt. implant



Cod	244-425-13
	244-425-13T
Deschidere maxima	7,5 mm
Fora inchidere	200 gf


 Titan  Otel pt. implant

*pentru inchidere permanenta



Cod	244-430-07	244-430-09	244-430-12	244-431-05	244-431-07	244-431-11
	244-430-07T	244-430-09T	244-430-12T	244-431-05T	244-431-07T	244-431-11T
Deschidere maxima	7,4 mm	8,4 mm	9,6 mm	5,8 mm	6,5 mm	7,2 mm
Forta inchidere	150 gf	150 gf	180 gf	150 gf	150 gf	180 gf

 Titan  Otel pt. implant




Cod	244-433-05	244-433-07	244-433-11
	244-433-05T	244-433-07T	244-433-11T
Deschidere maxima	4,9 mm	4,9 mm	4,9 mm
Forta inchidere	150 gf	150 gf	180 gf

 Titan  Otel pt. implant

CLIP ANEURISM YASARGIL FENESTRAT, 5 mm




*pentru inchidere permanenta




Cod	244-436-07	244-436-09	244-436-12	244-437-05	244-437-07	244-437-11
	244-436-07T	244-436-09T	244-436-12T	244-437-05T	244-437-07T	244-437-11T
Deschidere maxima	7,9 mm	9,1 mm	10,3 mm	6,5 mm	7,2 mm	7,8 mm
Forta inchidere	150 gf	150 gf	180 gf	150 gf	150 gf	180 gf

Titan Otel pt. implant



Cod	244-439-05	244-439-07	244-439-11	244-440-05	244-440-07	244-440-11
	244-439-05T	244-439-07T	244-439-11T	244-440-05T	244-440-07T	244-440-11T
Deschidere maxima	5,5 mm	5,5 mm	5,5 mm	7,9 mm	8,8 mm	10,0 mm
Forta inchidere	150 gf	150 gf	180 gf	150 gf	150 gf	180 gf


Titan Otel pt. implant



Cod	244-441-05	244-441-07	244-441-10	244-443-05	244-443-07	244-443-10
	244-441-05T	244-441-07T	244-441-10T	244-443-05T	244-443-07T	244-443-10T
Deschidere maxima	5,8 mm	5,8 mm	5,8 mm	5,8 mm	5,8 mm	5,8 mm
Forta inchidere	150 gf	150 gf	150 gf	150 gf	150 gf	150 gf


Titan Otel pt. implant

*pentru inchidere permanenta



Cod	244-445-05	244-445-07	244-445-10	244-447-05	244-447-07	244-449-05
	244-445-05T	244-445-07T	244-445-10T	244-447-05T	244-447-07T	244-449-05T
Deschidere maxima	7,0 mm	7,8 mm	8,8 mm	5,8 mm	6,2 mm	5,8 mm
Forta inchidere	150 gf	150 gf	150 gf	150 gf	150 gf	150 gf

 Titan Otel pt. implant




Cod	244-449-07	244-451-05	244-451-07	244-453-05	244-453-07
	244-449-07T	244-451-05T	244-451-07T	244-453-05T	244-453-07T
Deschidere maxima	6,2 mm	5,6 mm	5,6 mm	5,6 mm	5,6 mm
Forta inchidere	150 gf	150 gf	150 gf	150 gf	150 gf

 Titan Otel pt. implant

CLIP ANEURISM YASARGIL STANDARD




*pentru inchidere temporara




Cod	244-500-07	244-500-09	244-500-11	244-500-15	244-500-17	244-500-20
	244-500-07T	244-500-09T	244-500-11T	244-500-15T	244-500-17T	244-500-20T
Deschidere maxima	6,2 mm	7,0 mm	7,8 mm	9,2 mm	10,6 mm	11,4 mm
Forta inchidere	110 gf	90 gf	90 gf	90 gf	90 gf	110 gf

Titan Otel pt. implant



Cod	244-500-25	244-501-07	244-501-09	244-501-11	244-501-15	244-501-20
	244-500-25T	244-501-07T	244-501-09T	244-501-11T	244-501-15T	244-501-20T
Deschidere maxima	13,3 mm	6,0 mm	6,8 mm	7,5 mm	8,7 mm	11,0 mm
Forta inchidere	110 gf	110 gf	90 gf	90 gf	110 gf	90 gf


Titan Otel pt. implant



Cod	244-503-07	244-503-09	244-503-11	244-503-15	244-503-17
	244-503-07T	244-503-09T	244-503-11T	244-503-15T	244-503-17T
Deschidere maxima	5,5 mm	6,1 mm	6,5 mm	7,4 mm	9,9 mm
Forta inchidere	110 gf	110 gf	110 gf	110 gf	90 gf


Titan Otel pt. implant

*pentru inchidere temporara




Cod	244-505-07	244-505-09	244-505-11	244-510-05	244-510-07	244-510-09
	244-505-07T	244-505-09T	244-505-11T	244-510-05T	244-510-07T	244-510-09T
Deschidere maxima	5,4 mm	5,8 mm	6,2 mm	8,0 mm	7,9 mm	8,7 mm
Forta inchidere	110 gf	110 gf	110 gf	110 gf	90 gf	90 gf

Titan Otel pt. implant



Cod	244-510-12	244-510-17	244-510-20	244-511-05	244-511-07	244-511-09
	244-510-12T	244-510-17T	244-510-20T	244-511-05T	244-511-07T	244-511-09T
Deschidere maxima	10,2 mm	11,9 mm	13,0 mm	6,2 mm	6,5 mm	7,8 mm
Forta inchidere	90 gf	90 gf	90 gf	110 gf	110 gf	110 gf

Titan Otel pt. implant




Cod	244-511-11	244-513-05	244-513-07	244-513-10	244-515-08	244-515-10
	244-511-11T	244-513-05T	244-513-07T	244-513-10T	244-515-08T	244-515-10T
Deschidere maxima	8,4 mm	5,6 mm	5,6 mm	5,6 mm	5,4 mm	5,6 mm
Forta inchidere	110 gf	130 gf	130 gf	130 gf	130 gf	130 gf

Titan Otel pt. implant


CLIP ANEURISM YASARGIL STANDARD

*pentru inchidere temporara



Cod	244-517-05	244-519-09	244-519-11	244-521-07	244-521-11	244-523-09
	244-517-05T	244-519-09T	244-519-11T	244-521-07T	244-521-11T	244-523-09T
Deschidere maxima	7,2 mm	7,2 mm	8,2 mm	5,7 mm	8,0 mm	6,2 mm
Forta inchidere	90 gf	90 gf	90 gf	90 gf	90 gf	90 gf


 Titan Otel pt. implant



Cod	244-525-13
	244-525-13T
Deschidere maxima	7,5 mm
Forta inchidere	90 gf

 Titan Otel pt. implant


* pentru inchidere temporara



Cod	244-530-07	244-530-09	244-530-12	244-531-05	244-531-07	244-531-11
	244-530-07T	244-530-09T	244-530-12T	244-531-05T	244-531-07T	244-531-11T
Deschidere maxima	7,4 mm	8,4 mm	9,6 mm	5,8 mm	6,5 mm	7,2 mm
Forta inchidere	110 gf	110 gf	90 gf	110 gf	110 gf	110 gf

Titan

Otel pt. implant




Cod	244-533-05	244-533-07	244-533-11
	244-533-05T	244-533-07T	244-533-11T
Deschidere maxima	4,9 mm	4,9 mm	4,9 mm
Forta inchidere	110 gf	110 gf	110 gf

Titan

Otel pt. implant

**CLIP ANEURISM YASARGIL STANDARD,
PERFORATIE CIRCULARA 5,0 mm**


* pentru inchidere temporara



Cod	244-536-07	244-536-09	244-536-12	244-537-05	244-537-07	244-537-11
	244-536-07T	244-536-07T	244-536-12T	244-537-05T	244-537-07T	244-537-11T
Deschidere maxima	7,9 mm	9,1 mm	10,3 mm	6,5 mm	7,2 mm	7,8 mm
Forta inchidere	110 gf	100 gf	90 gf	110 gf	110 gf	110 gf

Titan


Otel pt. implant



Cod	244-539-05	244-539-07	244-539-11	244-540-05	244-540-07	244-540-11
	244-539-05T	244-539-07T	244-539-11T	244-540-05T	244-540-07T	244-540-11T
Deschidere maxima	5,5 mm	5,5 mm	5,5 mm	7,9 mm	8,8 mm	10,0 mm
Forta inchidere	110 gf	110 gf	90 gf	90 gf	90 gf	90 gf

Titan

Otel pt. implant




Cod	244-541-05	244-541-07	244-541-10	244-543-05	244-543-07	244-543-10
	244-541-05T	244-541-07T	244-541-10T	244-543-05T	244-543-07T	244-543-10T
Deschidere maxima	5,8 mm	5,8 mm	5,8 mm	5,8 mm	5,8 mm	5,8 mm
Forta inchidere	110 gf	110 gf	110 gf	110 gf	110 gf	110 gf

Titan

Otel pt. implant


* pentru inchidere temporara



Cod	244-545-05	244-545-07	244-545-10	244-547-05	244-547-07	244-549-05
	244-545-05T	244-545-07T	244-545-10T	244-547-05T	244-547-07T	244-549-05T
Deschidere maxima	7,0 mm	7,8 mm	8,8 mm	5,8 mm	6,2 mm	5,8 mm
Forta inchidere	110 gf	110 gf	110 gf	110 gf	110 gf	110 gf

Titan

Otel pt. implant



Cod	244-549-07	244-551-05	244-551-07	244-553-05	244-553-07
	244-549-07T	244-551-05T	244-551-07T	244-553-05T	244-553-07T
Deschidere maxima	6,2 mm	5,6 mm	5,6 mm	5,6 mm	5,6 mm
Forta inchidere	110 gf	110 gf	110 gf	110 gf	110 gf

Titan

Otel pt. implant

* pentru inchidere permanenta




Cod	244-460-03	244-460-05	244-462-03	244-462-05	244-462-07	244-463-03
	244-460-03T	244-460-05T	244-462-03T	244-462-05T	244-462-07T	244-463-03T
Deschidere maxima	3,3 mm	4,0 mm	3,3 mm	4,0 mm	4,6 mm	3,2 mm
Forta inchidere	110 gf	110 gf	110 gf	110 gf	110 gf	110 gf

Titan Otel pt. implant




Cod	244-463-04	244-465-04	244-465-05	244-467-04	244-467-05	244-467-07
	244-463-04T	244-465-04T	244-465-05T	244-467-04T	244-467-05T	244-467-07T
Deschidere maxima	3,8 mm	3,6 mm	3,8 mm	3,6 mm	4,0 mm	4,4 mm
Forta inchidere	110 gf	110 gf	110 gf	110 gf	110 gf	110 gf

Titan Otel pt. implant



Cod	244-469-04	244-471-04	244-471-05	244-473-05	244-475-06	244-476-04
	244-469-04T	244-471-04T	244-471-05T	244-473-05T	244-475-06T	244-476-04T
Deschidere maxima	3,6 mm	3,5 mm	4,0 mm	3,5 mm	6,0 mm	7,0 mm
Forta inchidere	110 gf	110 gf	110 gf	130 gf	130 gf	130 gf

Titan
Otel pt.
implant



Cod	244-476-07	244-477-07
	244-476-07T	244-477-07T
Deschidere maxima	5,7 mm	4,5 mm
Forta inchidere	130 gf	130 gf

Titan Otel pt. implant

*pentru inchidere permanenta



Cod	244-560-03	244-560-05	244-562-03	244-562-05	244-562-07	244-563-03
	244-560-03T	244-560-05T	244-562-03T	244-562-05T	244-562-07T	244-563-03T
Deschidere maxima	3,3 mm	4,0 mm	3,3 mm	4,0 mm	4,6 mm	3,2 mm
Forta inchidere	90 gf	70 gf	90 gf	70 gf	70 gf	90 gf

Titan


Otel pt. implant



Cod	244-563-04	244-565-04	244-565-05	244-567-04	244-567-05	244-567-07
	244-563-04T	244-565-04T	244-565-05T	244-567-04T	244-567-05T	244-567-07T
Deschidere maxima	3,8 mm	3,6 mm	3,8 mm	3,6 mm	4,0 mm	4,4 mm
Forta inchidere	70 gf	80 gf	70 gf	80 gf	70 gf	70 gf

Titan


Otel pt. implant



Cod	244-569-04	244-571-04	244-571-05	244-573-05	244-575-06	244-576-04
	244-569-04T	244-571-04T	244-571-05T	244-573-05T	244-575-06T	244-576-04T
Deschidere maxima	3,6 mm	3,5 mm	4,0 mm	3,5 mm	6,0 mm	7,0 mm
Forta inchidere	80 gf	80 gf	70 gf	70 gf	70 gf	70 gf

Titan

Otel pt. implant



Cod	244-576-07	244-577-07
	244-576-07T	244-577-07T
Deschidere maxima	5,7 mm	4,5 mm
Forta inchidere	70 gf	70 gf

Titan

Otel pt. implant